
EXPERIMENT 20

PREPARATION OF NYLON 66 - A CONDENSATION POLYMER

Structure

- 20.1 Introduction
 - Objectives
- 20.2 Polymers : an introduction
- 20.3 Principle
- 20.4 Requirements
- 20.5 Procedure
- 20.6 Result

20.1 INTRODUCTION

In the previous experiment you have learnt about preparation and use of azo dyes. In this process you have learnt about diazotisation reaction. In the present experiment you are going to prepare a sample of nylon 66 as an example of yet another important class of organic compounds viz. polymers.

In past few decades the developments in the field of polymers have dramatically changed our way of living. The way we dress ourselves, the way we pack our things, the way we decorate/furnish our household, everything has changed. Owing to their wide variety and versatility they have entered into practically every sphere of our day to day life. Polythene, used in carry bags; polyesters, the synthetic fabrics; poly vinyl chloride (PVC) used in water proof coverings; melamines, the unbreakable crockery; plastics, used for a whole range of materials have become part of our day to day vocabulary. The large number of available monomers and the ways in which they can combine together to form polymers provides for such a wide range of polymeric molecules,

In the preparation of nylon 66 you will learn about how small molecules "condense" together to make larger molecules i.e. condensation polymerisation. In the next section we have included a brief introduction to polymers, in general and nylon, in particular. In the next experiment, you will learn about another important class of organic compounds namely cosmetics and you will be preparing a sample of a face cream.


Objectives :

After performing these experiments you should be able to :

- prepare a sample of nylon 66
- explain condensation polymerisation, and
- list various applications of the polymer prepared.

20.2 POLYMERS : AN INTRODUCTION

The term 'POLYMERS' has been derived by the synthesis of two Greek words, "poly" meaning "many" and "meros" meaning "part". i.e. a molecule made up of many



(Nylong 66.)

The condensation of monomers to form the polymer can be done in a number of ways. You would be using the method of interfacial condensation. In this technique the polymerization is made to take place at the interface of an organic and an aqueous medium. This technique is good for the reactants with quite reactive functional groups, which is true in our case.

The diamine is dissolved in water and the diacid chloride is dissolved in an organic solvent like chloroform or carbon tetrachloride. The aqueous layer is carefully poured on the top of organic layer. Diamine being soluble in both the solvents diffuses from aqueous to organic layer and react with the reactive diacid chloride at interface and produces the insoluble polymer. This polymeric layer checks further diffusion of the diamine and the reaction stops. The by-product, HCl, formed during condensation diffuses back into the aqueous phase and gets neutralized by the base. The polymer can be removed with the help of a forcep or a glass rod and the process continues. The polymer is removed continuously till the reactant are exhausted.

At industrial scale production the polymer is not pulled out, as stated above, but the two solution are thoroughly agitated so as to form an emulsion. This leads to large increase in interfacial surface area. As a consequence the rate of polymerization increases and a good amount of polymer is obtained in short time.

20.4 REQUIREMENTS :

apparatus		chemicals
Round bottom flask 50 cm ³	1	Adipic acid
Air condenser	1	Hexamethylene diamine
Beaker 100 cm ³	2	Thionyl chloride
Beaker 400 cm ³	1	Dimethyl formamide
		Carbon tetrachloride
		Sodium hydroxide
		Alcohol


20.5 PROCEDURE :

Preparation of adipoyl chloride :

1. Set up the apparatus as shown in the figure. For this take a 50 cm³ round bottom flask and fit an air condenser on to it. Take a glass funnel and inset its stem into a piece of rubber tubing. In set a glass tube in the other end of the rubber tube and fit it on the air condenser as shown in the Fig. Given in margin. Put the inverted funnel into a beaker containing water or NaOH solution.

2. Take 1g of adipic acid in 3-4 drops of dimethyl formamide in the flask and add 1 cm³ of thionyl chloride dropwise with constant stirring.

3. Heat the flask on water bath for about 15 minutes. By this the evolution of gas ceased and the solid will disappear.



**Preparation of Nylon 66 - a
Condensation Polymer**

4. Allow the flask to cool a little and then add about 30cm³ of CCl₄ to it. Mix thoroughly to dissolve the product obtained.
5. Transfer this solution to a 100 cm³ beaker. Wash the flask thoroughly with additional 20 cm³ of CCl₄ and transfer the washings to the beaker.
6. In a separate 100 cm³ beaker take 1.1 g of hexamethylene diamine and 0.75 g (6-7 pallets) of NaOH in 25 cm³ of water and mix to dissolve them.
7. Carefully transfer the aqueous solution of hexa methylene diamine to the beaker containing adipoyl chloride.
8. You will observe the formation of a film of nylon 6, 6, at the interface of the two liquids. Carefully insert a glass rod or a copper wire into the solution and pull out the polymer formed. You may even use a forcecup to do so.
9. Wrap the polymer around a clean test tube as shown in the diagram. Rotate the test tube to pull more and more of nylon.

CAUTION ! Avoid handling the polymer with hands at this stage because the reactants specially adipoyl chloride and the solvent (CCl₄) both are harmful.

Do not pour the unreacted polymerization mixture in the sink. Stir it with a glass rod till there is no more polymerization. You may recover this crop of the polymer also and discard the rest of the solvent as usual.

10. Wash the polymer thoroughly with water. For this place the test tube, with polymer wrapped around, under the running tap water for about 5 minutes. Alternatively you may first wash the polymer with 50 % aqueous alcohol (alcohol:H₂O ::1:1) followed by tap water.
11. Dry the polymer in air or in the folds of filter paper. Weigh it and report the yield.
12. Submit a sample of the polymer so prepared to your counsellor.

20.6 RESULT

..... g of nylon 66 was obtained from g of adipic acid.

Put 1-2 crystals of azobenzene, if available, to this solution and mix to dissolve. This will enhance the visibility of the interface of two liquids.

OR

You may add a drop of phenolphthalein to the aqueous layer containing hexamethylenediamine.

You should add the aqueous solution slowly along the walls of the beaker. Taking care, not to allow mixing of the solvents.